


Kanak Raj, Anurag P. Srivastava, Thummala Praveen, Rajiv Bharadwaj, Sunit Maity, Vaishali Sinha and Sridevi Khambhampaty, Biopharmaceutical Development Department, Syngene International Ltd, Biocon Park, Bengaluru Karnataka – 560099, E. Mail: kanak.raj@syngeneintl.com

### Introduction


Plasmid DNA is a new generation biotechnology product (gene medicines and DNA vaccines) that is just beginning to enter the marketplace. Plasmid DNA vectors may find application as preventive or therapeutic DNA vaccines for viral, bacterial, parasitic diseases or for other indications such as cancer and gene therapy products [1]. However, large scale production of pDNA is challenging, due to effect of production parameters on the extent of plasmid supercoiling, generational loss and instability [2]. Syngene scientists were able to overcome these challenges with their extensive experience and novel ideas, resulting in a pDNA production platform process with high pDNA yield and purity. Microbial Upstream team at Syngene has developed a high cell density fed-batch fermentation platform process to produce pDNA at high titers (>1 g/L), while downstream scientists have overcome the challenges by implementing a custom-made cell lysis vessel, which allowed to control critical parameters such as lysis pH, agitation, and temperature. The Process Development was supported by state-of-the-art analytical capability at Syngene.

### Process


#### Overview




#### Downstream



#### Upstream




#### Analytical Methods For pDNA



### Results

#### Upstream



#### Downstream

#### DSP FLOW DIAGRAM

```

graph TD
    A[5% Cell suspension in solution A] --> B[Equal volume of soln. B followed by 5 min incubation @8-10°C]
    B --> C[Added soln. C (30min @8-10°C)]
    C --> D[Centrifugation]
    D --> E[CaCl2 precipitation]
    E --> F[UFDf]
    F --> G[Anion Exchange Chromatography]
    G --> H[Affinity Chromatography]
    H --> I[Gel filtration Chromatography]
        
```

#### AUTOMATED LYSIS SYSTEM (ALS)

**Benefits of ALS:**

- Automated lysis system helps in performing continuous lysis with less shear under controlled condition of critical parameters during lysis.
- ALS provides good yield of supercoiled pDNA.

#### AGAROSE GEL ELECTROPHORESIS

Attributes	Percentage
Recovery	> 25%
Purity	> 90%

**Plasmid DNA Downstream Platform Process:**

- Syngene has developed a scalable platform process for the plasmid DNA production.
- The developed process worked equally well with different plasmid sizes (2-14 kb).
- Process resulted in yield of 25-40% enriched with the supercoiled form of plasmid DNA

### Conclusion

After successful planning and experimentation, Syngene scientists were able to develop a platform process for plasmid DNA production, with high titers > 1 g/L at Upstream and >25% recovery at Downstream. The platform process has been successfully tested across 2L, 10L and 50L production scales, for both high and low copy number plasmids. Syngene pDNA platform process has served and continues to serve the needs of a diversity of Syngene clients.

### Acknowledgement

The authors would like to thank Syngene International Ltd., for funding this project

### References

- Mark A. Eiteman and Elliot Altman, Overcoming acetate in *Escherichia coli* recombinant protein fermentations, 2006, TRENDS in Biotechnology Vol.24 No.11
- James A. Williams, Aaron E. Carnes, Clague P. Hodgson, Plasmid DNA vaccine vector design - Impact on efficacy, safety and upstream production, 2009, Biotechnology Advances, 27, 353-370